Sino-lipo2000 Transfection Reagent

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2. Transfect cells according to the following chart. Volumes are given on a per-well basis. Each reaction mix is sufficient for triplicate (96-well), duplicate (24-well), and single well (6-well) transfections, and accounts for pipetting variations. Adjust the amounts of components according to your tissue culture format.

Timeline		Steps	Procedure Details			
Day 0	1	Seed cells to be 70~90% confluent at transfection	Component	96-well	24-well	6-well
	2	Dilute four amounts of Lipo2000 Reagent in Opti-MEM® Medium	Adherent cells	$1-4 \times 10^4$	$0.5-2 \times 10^5$	$0.25-1 \times 10^6$
	3	Dilute DNA in Opti-MEM® Medium	Opti-MEM® Medium	25 $\mu\text{L} \times 4$	50 $\mu\text{L} \times 4$	150 $\mu\text{L} \times 4$
Day 1	4	Add diluted DNA to diluted Lipo 2000 Reagent (1:1 ratio)	Lipo2000 Reagent	① 1, 1.5, 2, 2.5 μL	② 2, 3, 4, 5 μL	③ 6, 9, 12, 15 μL
	5	Incubate	Opti-MEM® Medium	125 μL	250 μL	700 μL
	6	Add DNA-lipid complex to cells	DNA (0.5-5 $\mu\text{g}/\mu\text{L}$)	2.5 μg	5 μg	14 μg
	7	Visualize/analyze transfected cells	Diluted DNA Total	25 μL	50 μL	150 μL
			Diluted Lipo 2000 Reagent	25 μL	50 μL	150 μL
			Incubate for 5 minutes at room temperature.			
			Component	96-well	24-well	6-well
			DNA-lipid complex per well	10 μL	50 μL	250 μL
			Final DNA used per well	100 ng	500 ng	2500 ng
			Final Lipo 2000 Reagent used per well	0.2-0.5 μL	1.0-2.5 μL	5.0-12.5 μL
			Incubate cells for 1-3 days at 37°C. Then analyze transfected cells.			

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